

Product name	Microcentrifuge Tubes
Description / Specification	
Materials	Polypropylene (PP) Meets USP, Class VI standard
Temperature range	stored at room temperature
Shelf life	3 Years after date of production (ensure package is in good)
Sterilization	Yes Sterilized by Irradiation , Sterility Assurance Level: SAL=10 ⁻⁶ . The product has been irradiated and dosimetrically released based on ANSI/AAMI/ISO 11137
Clolour	Clear
Pyrogens	Non-Pyrogenic
RNase/DNase	DNase/RNase free
Endotoxin	endotoxin test<0.005EU/ML
Performance Testing	Each manufacturing lot is sampled and tested in accordance with standard operating procedures. Appearance inspection: qualified Packaging inspection: qualified Sealing test: qualified Freeze testing: qualified
Features	
Basic features	1. None of DNA\RNA enzyme, endotoxin<0.005EU/ML . 2. Using imported PP materials, resistant to chemical corrosion and low temperature, compliant to USP Class VI standard. 3. Special safelock cap design provides secure lid sealing. 4. Tube cap can be opened or closed with one hand, easy to operate.
Temperature range	-80°C-121°C
Centrifugation, max. RCF	20000xg
Volumes	0.6ml,1.5ml,2.0ml
Packaging	Bulk packed

Product Description			
Cat.No.	Product Description	Pack	Case
311001	0.6mL Microcentrifuge Tube, Conical, Lock Cap, Non-sterile	1000	10000
311002	1.5mL Microcentrifuge Tube, Conical, Lock Cap, Non-sterile	500	5000
311004	2.0mL Microcentrifuge Tube, Conical, Lock Cap, Non-sterile	500	5000

DISCLAIMER: The description of a certain product can only be considered as a guidance, because its performance ultimately depends on what the product is used for. Very often performance studies are indispensable.

Qingdao Haier Biomedical Co., Ltd
Add: No.280 Fengyuan Road, High-tech Zone, Qingdao P.R. China
Tel:(86)532-88935593
Website: www.haiermedical.com