

Product name	Filter Pipette Tips
Description / Specification	
Materials	Polypropylene (PP). Meets USP, Class VI standard
Temperature range	stored at room temperature
Shelf life	3 Years after date of production (ensure package is in good)
Sterilization	Yes (for Boxed tips) Sterilized by E-beam, Sterility Assurance Level: SAL=10 ⁻⁶ . The product has been irradiated and dosimetrically released based on ANSI/AAMI/ISO 11137
Pyrogens	Non-Pyrogenic
RNase/DNase	DNase/RNase free
Quality Control	Each manufacturing lot is sampled and tested in accordance with standard operating procedures. Raw Material-Control: physical testing Product-Control: testing of attributive and variable characteristics in accordance with the valid specification
Features	
Basic features	1. Both tips and boxes are autoclavable. 2. Filter Pipette Tips are filtered through porous tissue to ensure optimal performance and is extremely hydrophobic, forming a robust barrier against aerosols and eliminating the risk of cross-contamination of the sample. 3. Haier pipette tips are compatible with most of the popular

Qingdao Haier Biomedical Co., Ltd

Add: No.280 Fengyuan Road, High-tech Zone, Qingdao P.R. China

Tel:(86)532-88935593

Website: www.haiermedical.com

	pipettes on the market. 4. Ultra-clear polypropylene good for observation.
Volumes	10/20/50/100/200/300/1000/1250 μ L
Packaging	Bagged and racked

Product Description							
Cat.No	Volume(μ L)	Type	Sterile	Colour	Pack	Low retention	Case
112005	10	Bagged	×	Clear	1000	×	10000
111006	10	Racked	√	Clear	96	×	4800
111007	10	Racked	√	Clear	96	√	4800
112048	10(Extra Long)	Bagged	×	Clear	1000	×	10000
111049	10(Extra Long)	Racked	√	Clear	96	×	4800
111050	10(Extra Long)	Racked	√	Clear	96	√	4800
112008	20	Bagged	×	Clear	1000	×	10000
111009	20	Racked	√	Clear	96	×	4800
111010	20	Racked	√	Clear	96	√	4800
112011	50	Bagged	×	Clear	1000	×	10000
111012	50	Racked	√	Yellow	96	×	4800
111013	50	Racked	√	Clear	96	×	4800

Qingdao Haier Biomedical Co., Ltd

Add: No.280 Fengyuan Road, High-tech Zone, Qingdao P.R. China

Tel:(86)532-88935593

Website: www.haiermedical.com

111014	50	Racked	√	Clear	96	√	4800
112015	100	Bagged	×	Clear	1000	×	1000 0
111016	100	Racked	√	Clear	96	×	4800
111017	100	Racked	√	Clear	96	√	4800
112022	200	Bagged	×	Clear	1000	×	1000 0
111023	200	Racked	√	Clear	96	√	4800
111024	200	Racked	√	Clear	96	√	4800
111031	200(Extra Long 59mm)	Racked	√	Clear	96	×	4800
122033	200(Extra Long 59mm)	Racked	√	Clear	96	√	4800
111035	200(Extra Long 78mm)	Racked	√	Clear	96	×	4800

Qingdao Haier Biomedical Co., Ltd

Add: No.280 Fengyuan Road, High-tech Zone, Qingdao P.R. China

Tel:(86)532-88935593

Website: www.haiermedical.com

111037	200(Extra Long 78mm)	Racked	√	Clear	96	√	4800
112041	300(Extra Long)	Bagged	×	Clear	1000	×	5000
111042	300(Extra Long)	Racked	√	Clear	96	×	4800
111043	300(Extra Long)	Racked	√	Clear	96	√	4800
112057	1000	Bagged	×	Blue	1000	×	5000
112058	1000	Bagged	×	Clear	1000	×	5000
111059	1000	Racked	√	Blue	100	×	5000
111060	1000	Racked	√	Clear	100	√	5000
112066	1250	Bagged	×	Blue	1000	×	5000
112067	1250	Bagged	×	Clear	1000	×	5000
111068	1250	Racked	√	Blue	96	×	4800
111069	1250	Racked	√	Clear	96	×	4800
111070	1250	Racked	√	Clear	96	√	4800

DISCLAIMER: The description of a certain product can only be considered as a guidance, because its performance ultimately depends on what the product is used for. Very often performance studies are indispensable.

Qingdao Haier Biomedical Co., Ltd

Add: No.280 Fengyuan Road, High-tech Zone, Qingdao P.R. China

Tel:(86)532-88935593

Website: www.haiermedical.com