

Product name	Centrifuge Tubes
Description / Specification	
Materials	Tube: Polypropylene (PP) Cap: High Density Polyethylene (HDPE)
Temperature range	stored at room temperature
Shelf life	3 Years after date of production (ensure package is in good)
Sterilization	Yes Sterilized by E-beam, Sterility Assurance Level: SAL=10 ⁻⁶ . The product has been irradiated and dosimetrically released based on ANSI/AAMI/ISO 11137
Pyrogens	Non-Pyrogenic
RNase/DNase	DNase/RNase free
Colour	Tube: clear Cap: blue
Performance Testing	Each manufacturing lot is sampled and tested in accordance with standard operating procedures. Appearance inspection: qualified Sealing test: qualified Packaging inspection: qualified
Features	
Basic features	Non-Pyrogenic, DNase/Rnase free
Temperature range	-80℃ -121℃
Volumes	15/50 mL
Centrifugation, max. RCF	16000 xg
Packaging	Bulk packed in zip-Lock bag

Product Description				
Cat.No.	Volume(mL)	Max.RCF (xg)	Pack	Case
311005	15	16000	25	500
311006	50	16000	25	500

DISCLAIMER: The description of a certain product can only be considered as a guidance, because its performance ultimately depends on what the product is used for. Very often performance studies are indispensable.