

Product name	Serological Pipette
Description / Specification	
Materials	100% Virgin Polystyrene (PS).
Temperature range	stored at room temperature
Shelf life	3 Years after date of production
Sterilization	Yes Sterilized by E-beam, Sterility Assurance Level: SAL=10 ⁻⁶ .The product has been irradiated and dosimetrically released based on ANSI/AAMI/ISO 11137
Pyrogens	Non-Pyrogenic
RNase/DNase	DNase/RNase free
Performance Testing	Each manufacturing lot is sampled and tested in accordance with standard operating procedures. Appearance inspection:qualified Stress test:qualified Sealing test: qualified Packaging inspection: qualified
Features	
Basic features	Edged color-coded band for identifying Pipettes of different volumes. Negative graduations for extra volume. Polyolefin fiber filter to reduce contamination.
Temperature range	Approx. 0℃to 40℃
Volumes	1 mL/ 2 mL /5 mL /10 mL/ 25 mL /50 mL /100mL
Packaging	Individually wrapped in paper-plastic bags

Product Range				
Cat.No.	Volume(mL)	Color Ring	Pack	Case
692209	1	Yellow	1	1000
692210	2	Green	1	1000
692211	5	Blue	1	200
692212	10	Orange	1	200
692213	25	Red	1	250
692214	50	Purple	1	100
692215	100	Black	1	120

DISCLAIMER: The description of a certain product can only be considered as a guidance, because its performance ultimately depends on what the product is used for. Very often performance studies are indispensable.